

LIST OF KNOWN COMMERCIALY-AVAILABLE ANTIGEN-DETECTING MALARIA RDTs
Information for national public health services and UN Agencies wishing to procure RDTs
 (This list is not based on any criteria of product quality)

See notes on purchase and use of RDTs to guide use of this list, <http://www.wpro.who.int/rdt/link4.cfm>

Disclaimer: Inclusion in this list does not in any way imply an endorsement, certification, warranty of fitness or recommendation by WHO of any company or product for any purpose, and does not imply preference over products of a similar nature that are not mentioned. WHO does not furthermore warrant that: (1) the list is complete and/or that the information provided is or will be always accurate; and/or that (2) the products listed are of acceptable quality, have obtained regulatory approval in any country, or that their use is otherwise in accordance with the national laws and regulations of any country, including but not limited to patent laws. The list is based on evidence of quality manufacturing, through provision of independent certification of compliance with ISO13485:2003 (or US FDA 21 CFR 820) and a stability testing protocol provided by the companies to WHO, and whilst requesting this evidence from manufacturers as a prerequisite for listing of products. WHO can not guarantee the accuracy of this information. WHO disclaims any and all liability and responsibility for any injury, death, loss, damage, or other prejudice of any kind whatsoever that may arise as a result of, or in connection with the procurement, distribution and use of any product included on the list. WHO recommends that the products should only be used by individuals trained

Products	Product name	Plasmodium species targeted (F = P falciparum V = P vivax O = P ovale M = P malariae P = pan-specific)	Target Antigen	Format	Packaging	All materials can be included? (Lancets, swabs, wells if required)	Shelf Life (months)	Storage temperature and stability data on which Shelf Life is based (Acc = Accelerated or RT = Real Time)	Contact	Address & Contact Nos.	Email	URL
Access Bio	CareStart Malaria pLDH PAN Test	P	pLDH	Cassette	Individual or Bulk	Yes	18	4-30 C, RT	Jaeun Jung	2033 Route 130 Unit-H Mormouth Junction, NJ 08852 USA Tel: +1 732-297-2222 x102 Fax: +1 732-297-3001	jjung@accessbio.net	www.accessbio.net
	CareStart Malaria pLDH 3-Line Test	F, P, V	pLDH	Cassette	Individual or Bulk	Yes	18	4-30 C, RT				
	CareStart Malaria pLDH/HRP2 Com	F, P	pLDH/HRP2	Cassette	Individual or Bulk	Yes	18	4-30 C, Acc				
	CareStart Malaria HRP2 Test	F	HRP2	Cassette	Individual or Bulk	Yes	18	4-30 C, Acc				
	CareStart Malaria pLDH Pf Test	F	pLDH	Cassette	Individual or Bulk	Yes	18	4-30 C, RT				
ACON Laboratories, Inc.	Malaria Pf Rapid Test Device	F	HRP2	Cassette	1 testbox, Individual	Test device, Buffer, Swab, lancet, dropper, packaging insert	24	2-30 C(Acc data)	Mr Lee Mariano Business Development Director	4108 Sorrento Valley Blvd., San Diego, CA 92121, U.S.A. Tel. no. 1-858-875-8000 (tel) Fax no. 1-858-535-2035	mariano@aconlabs.com info@aconlabs.com	www.aconlabs.com
	Malaria Pf Rapid Test Device	F	HRP2	Cassette	40 tests/box, bulk	Test device, Buffer, Swab, dropper, packaging insert	24					
AZOG	AZOG Malaria Pf (HRP2)/Pv (Pan) Antigen detection Test Device	F, P	HRP2, pLDH	Cassette	25tests/box, individual or bulk	Yes (5ul calibrated pipette, retractable lancet,alcohol)	18	RT (4-30)	Dr. Azubuke Ogata	1011 US Highway 22 West, Phillipsburg, NJ 08865 USA Tel: 908-213-2900	info@azoglabs.com	www.azoglabs.com
All Diag	PALUTOP	F	HRP2 pLDH specific for P. vivax pLDH common for 4 species	Cassette	Individual	Yes	18	4-30 C, RT	Dr Thomas Lamy Chief Executive Officer	12 rue Etienne Bugault, BP6 67038 Strasbourg CEDEX 2, France Tel: +33-3-8878-8088 Fax: +33-3-8878-7678	info@alldiag.com	www.alldiag.com
	PALUTOP +4	F, V, P	HRP2	Cassette	Individual	Yes	18	4-30 C, RT				
Ameritek, U.S.A. (Ameritech Diagnostic Reagent (JiaXing) Co., Ltd)	sBest One Step Rapid Malaria Test	F	HRP2	Dipstick	50 test/ pack Individual and Bulk	Yes	30	4-28 C, RT	K.C. Yee, M. D., Ph. D President and CEO Brad Butler Account Manager	AmeriTek USA 125 130th St., SE Suite #200 Everett, WA 98208 Tel: +1-425-379-2590 Fax: no. 425-379-2624	kcyee@ameritek.org info@ameritek.org www.ameritek.org	
	sBest One Step Malaria Test	F	HRP2	Cassette	30 test/pack Individual and Bulk	Yes	30	4-28 C, RT				
	sBest One Step Malaria Test	V	pLDH	Dipstick	30 test/pack Individual and Bulk	Yes	30	4-28 C, RT				
	sBest One Step Malaria Test	V	pLDH	Cassette	30 test/pack Individual and Bulk	Yes	30	4-28 C, RT				
Amgenix International, Inc.	OnSight - Malaria Pf (Test Device)	F	HRP2	Cassette	1.5, 10, 25, 50, 100 tests/box.	Yes	24	4-30 C, RT	A. Joan Director	3444 Pinotn Ct., San Jose, CA 95148 United States of America Tel. no. (408) 270-6241 Fax no. (408) 270-6214	amgenix@aol.com ajjoan@amgenix.com	www.amgenix.com
	OnSight - Malaria Pf (Test Strip)	F	HRP2	Dipstick	1.5, 10, 25, 50, 100 tests/box.	Yes	24	4-30 C, RT				
	OnSight - ParaQuick (Pan, Pf)	P, F	pLDH/HRP2	Cassette	1.5, 10, 25, 50, 100 tests/box.	Yes	24	4-30 C, RT				
	OnSight - ParaQuick-2 (Pv, Pf)	V, F	pLDH/HRP2	Cassette	1.5, 10, 25, 50, 100 tests/box.	Yes	24	4-30 C, RT				
	OnSight - ParaQuick-3 (Pan, Pv, Pf)	P, V, F	pLDH/HRP2	Cassette	1.5, 10, 25, 50, 100 tests/box.	Yes	24	4-30 C, RT				
	OnSight - PanScreen (Pan) (Test Device)	P	pLDH	Cassette	1.5, 10, 25, 50, 100 tests/box.	Yes	24	4-30 C, RT				
	OnSight - PanScreen (Pan) (Test Strip)	P	pLDH	Dipstick	1.5, 10, 25, 50, 100 tests/box.	Yes	24	4-30 C, RT				
Bhat Biotech	Maleriscan Malaria PTPV 3 line card test (antibody)	F, V, O, M	MSP, CSP	Cassette	Individual 50 tests/box	Yes	18	2-30 C, RT	Dr Shama Bhat Managing Director Tel: 91-80-2783 3620; 91-80-2783 3473 Fax: 91 80-2783 3621	No 11-A, 4th Cross, Veerasandra Indl. Area, Electronics City, Bangalore-560 100, Karnataka, India	bhatbiotech@vsnl.net	www.bhatbiotech.com
	Maleriscan Malaria PTPV 3 line antigen test	F, V, O, M	HRP, LDH	Cassette	Individual 50 tests/box	Yes	18	2-30 C, RT				
Binox, Inc.	NOW Malaria	P	HRP2 for Pf	Card	Individual	No	14	15-30 C, RT	Brian Conibere International Marketing Manager - Unipath US Contact: Jonathon Spears	Priority Business Park Bedford, MK44 3UP United Kingdom Tel: +1 207-772-3988 ext. 107 Fax: +1 207-871-5751	brian.conibere@unipath.com jonathon.spears@unimed.com	www.unipath.com
Biosynex	Immunoquick Malaria	F	HRP2	Dipstick	25 tests/box, Individual	Yes	18	4-30 C, Acc	Dr Thierry Paper General Manager	12 rue Etienne Bugault, BP6 67038 Strasbourg CEDEX 2, France	paper@biosynex.com	www.biosynex.com
	Immunoquick Malaria +4	F, P	HRP2/pLDH	Dipstick	25 tests/box, Individual	Yes	18	4-30 C, Acc				
Biotech Laboratories Ltd.	Malaria Pf rapid device	F	HRP2	Cassette	25 tests/pack, Individual	Yes	18	4-30 C, RT	Mike Ellis Sales Steve Birch Sales	32 Anson Road Marlsham Health Ipswich, Suffolk, IP5 3RG United Kingdom Tel: +44 1473 612168 Fax: +44 1473 611476	mike.ellis@biotec.com steve.birch@biotec.com	www.biotech.com
	Malaria Pv/Pf rapid device	V, F	pLDH/HRP2	Cassette	25 tests/pack, Individual	Yes	18	4-30 C, RT				
	Malaria Pan/Pf rapid device	P, F	HRP2	Cassette	25 tests/pack, Individual	Yes	18	4-30 C, RT				
	Malaria Pan/Pv/Pf rapid device	P, V, F	pLDH/HRP2	Cassette	25 tests/pack, Individual	Yes	18	4-30 C, RT				
Cellabs Pty. Ltd.	Rapimal Dipstick	F	HRP2	Dipstick	Individual	Yes	18	RT	Dr Hubert Mazure Marketing Manager	Unit 7, 27 Dale St. Brookvale, NSW 2100 Australia Tel: +61-2-9905-0133 Fax: +61-2-9905-6426	enquiries@cellabs.com.au	www.cellabs.com.au
	Rapimal Cassette	F	HRP2	Cassette	Individual	Yes	18	RT				
Chemoquip Ltd.	Eurocheck Pf	P	HRP-2	Cassette	Individual	Yes	24	4-40 C	Mr Ferroz Niwab	Ngara Road, P. O. Box 32565, 00600, Nairobi, Kenya Tel: 3749383/7, 3746511, 3747951 Fax: 3750274, 3746661	cnj@chemoquip.com	www.chemoquip.com
	Eurocheck Pf	P	HRP-2	Dipstick	Individual	Yes	24	4-40 C				

LIST OF KNOWN COMMERCIALY-AVAILABLE ANTIGEN-DETECTING MALARIA RDTs
Information for national public health services and UN Agencies wishing to procure RDTs
 (This list is not based on any criteria of product quality)

See notes on purchase and use of RDTs to guide use of this list, <http://www.wpro.who.int/rdt/link4.cfm>

Disclaimer: Inclusion in this list does not in any way imply an endorsement, certification, warranty of fitness or recommendation by WHO of any company or product for any purpose, and does not imply preference over products of a similar nature that are not mentioned. WHO does not furthermore warrant that: (1) the list is complete and/or that the information provided is or will be always accurate; and/or that (2) the products listed are of acceptable quality, have obtained regulatory approval in any country, or that their use is otherwise in accordance with the national laws and regulations of any country, including but not limited to patent laws. The list is based on evidence of quality manufacturing, through provision of evidence of independent certification of compliance with ISO13485:2003 (or US FDA 21 CFR 820) and a stability testing protocol provided by the companies to WHO, and whilst requesting this evidence from manufacturers as a prerequisite for listing of products. WHO can not guarantee the accuracy of this information. WHO disclaims any and all liability and responsibility for any injury, death, loss, damage, or other prejudice of any kind whatsoever that may arise as a result of, or in connection with the procurement, distribution and use of any product included on the list. WHO recommends that the products should only be used by individuals trained

Products	Product name	Plasmodium species targeted (F = P falciparum V = P vivax O = P ovale M = P malariae P = pan-specific)	Target Antigen	Format	Packaging	All materials can be included? (Lancets, swabs, wells if required)	Shelf Life (months)	Storage temperature and stability data on which Shelf Life is based (Acc = Accelerated or RT = Real Time)	Contact	Address & Contact Nos.	Email	URL
Core Diagnostics Ltd.	Core Malaria, cat: MAL-190020	F	HRP2	Cassette	25 tests/pack, Individual	Yes	22	4-30 C, RT	Michael Ittah	Aspect Court Road, 4 Temple Row, Birmingham B2 5HG United Kingdom Tel: +44-121-609-4720 Fax: +44-121-609-4721	sales@corediag.com	www.core-diagnostics.com
	Core Malaria, cat: MAL-190024	F, P	pLDH/HRP2	Cassette	25 tests/pack, Individual	Yes	22	4-30 C, RT				
	Core Malaria, cat: MAL-190022	V, F	pLDH	Cassette	25 tests/pack, Individual	Yes	22	4-30 C, RT				
	Core Malaria, cat: MAL-190026	F, V, P	pLDH/HRP2	Cassette	25 tests/pack, Individual	Yes	22	4-30 C, RT				
Diagnostic Automation (Cortez Diagnostics, Inc.)	Malaria Pf / Vivax	F, P - V, M, O		Cassette	Bulk	Yes	24	4-28 C	Marcia Muller Technical Support Manager	2361 Craftsman Road, Suite E/F, Calabasas, California 91302 USA Tel: +1 818-591-3030 Fax: +1 818-591-8383	technicalsupport@rapidtest.com	www.rapidtest.com
DIALAB	Malaria P. falciparum Cassette	F	HRP2	Cassette	Kitboxed Individually packed	Yes	24	2-30 C	Martin Herr	Z-NO Sud, Hondastrasse, Objekt M55, A-2351 Wr. Neudorf Tel: +43(0)2236 660910-0	herr@dalab.at	www.dialab.at
Diamed AG Diagnostic and Medical Products	OptiMal IT	F, P	pLDH	Cassette	Individual and Bulk	Yes	12	2-30 C, RT	Dr Patrick Jacquier Head of Parasitology and Infectious Diseases	Diamed SA, 1785, Cressier sur Morat, Switzerland Tel: +4126 674 5111 Fax: +4126 674 5145	p.jacquier@diamed.ch parasitology@diamed.ch	www.diamed.ch/optimal
	OptiMal 48	F, P	pLDH	Dipstick	Individual and Bulk	Yes	12	2-30 C, RT				
Fortress Diagnostics Ltd	PI000010	F	HRP2	Cassette	10/30 tests/pack	Yes	18	15-25 C, Acc	Dr Morteza Afraabi	Unit 2C Antrim Technology Park Belfast Road, Antrim BT41 1QS United Kingdom Tel: +44 2894 487676	info@fortressdiagnostics.com	www.fortressdiagnostics.com
	One Step Malaria Blood Test	F		Cassette	Individual/bulk	Yes. For individual pack, lancet, swab, buffer, dropper, test device and package insert are included. For bulk packaging, device, buffer and dropper are included if required.	24	4-30 C, Acc				
Guangzhou Wondfo Biotech Co. Ltd.	One Step Malaria Blood Test	F		Dipstick	Individual/bulk	Yes. For individual pack, lancet, swab, buffer, test device and package insert are included. For bulk packaging, device, buffer and dropper are included if required.	24	4-30 C, Acc	Ms Ivy Cheng Ms Sherry Liu	Wondfo Scientific Park South China University of Technology Guangzhou, P.R. China Tel: no. +86 20 705 7121 Fax no. +86 20 8711 1434	ivy@wondfo.com.cn sherry@wondfo.com.cn sales@wondfo.com.cn	www.wondfo.com.cn
	Hexagon Malaria	F	HRP2	Dipstick	20 tests/box, Individual	Lysing buffer, test tubes, capillary pipettes	18	2-30 C, RT	Dr Michael Dreher Technical Director	Max-Planck-Ring 21 65205 Wiesbaden, Germany Tel: +49 6122 9988 0 Fax: +49 6122 9988 100	human@human.de	www.human.de www.human-de.com
Hexagon Malaria	F	HRP2	Cassette	25 tests/box, Individual	Lysing buffer, capillary pipettes	18	2-30 C, RT					
Hexagon Malaria Combi	F, P	HRP2	Dipstick	20 tests/box, Individual	Lysing buffer, test tubes, capillary pipettes	18	2-30 C, Acc					
Hexagon Malaria Combi	F, P	HRP2	Cassette	25 tests/box, Individual	Lysing buffer, capillary pipettes	18	2-30 C, Acc					
ICT Diagnostics	ICT Malaria Pf	F	HRP2	Dipstick	25 tests/box	Yes	24	2-30 C, RT	Mr Russell Glanz Director	P.O. Box 912, Noordhoek 7985 Cape Town, South Africa Tel/Fax: +27-21-789-2979	russellag@icon.co.za	
	ICT Malaria Pf	F	HRP2	Cassette	25 or 100 tests/box	Yes	24	2-30 C, RT				
	ICT Malaria Combo	F, P, V	HRP2	Dipstick	25 tests/box	Yes	24	2-30 C, RT				
	ICT Malaria Combo	F, P, V	HRP2	Cassette	25 or 100 tests/box	Yes	24	2-30 C, RT				
IND Diagnostics	One Step Malaria Antigen Test	P	pLDH	Dipstick	25 tests/Kit, Individual or bulk	All materials can be included (no swab)	12	room temperature, Acc and RT data in progress	Jason Peng Shuli Wang	1629 Fosters Way, Delta B.C. V3M 6S7 Canada	ra@ind.ca Shuli@ind.ca	
Innovatek Medical Inc.	Quickstick Malaria Antigen Test	P,F	pLDH	Dipstick		Alcohol Swab, Lancets, Gauze or cotton balls not provided	2-30 C, Acc	Mr King F Hui	Innovatek Medical, 3-1600, Derwent Way, Delta B.C. Canada, V3M 6M5	khui@innovatekmed.com		
InTec PRODUCTS, INC.	ADVANCED QUALITY™ One Step Malaria (pL)	F	HRP2	Card/strip	40T/Pack Individual and Bulk	No	24	2-30 C	Mr Huang Qing Bo	InTec PRODUCTS, INC., 322, Xingjiang Road, Jinyang Industrial Area, Haining, 311022, Zhejiang, P.R. China	intecproducts@asintec.com	www.intec.com
	ADVANCED QUALITY™ MALARIA (pL) POCT	F	HRP2	Card	1T/Pack and Bulk	Test device, Buffer, Swab, lancet, dropper, lancets, sterile wipes or pads, timer not provided	24	2-30 C				
Inverness Medical Innovations	Sinax Now Malaria	F, P	HRP2/pan malaria antigen	Cassette			2-37 C	Ms Amy Yorston	Inverness Medical Innovations Australia Pty Ltd, 532 Seventeen Mile Rocks Road Beverly Hills QLD 4003	Amy_Yorston@panbio.com.au	www.invernessmedicalpd.com.au	
J. Mitra & Co. Pvt. Ltd.	Advantage P.I. Malaria Card	F	HRP2/	Cassette	25/50/100 tests/pack Individual	No	24	4-30 C, Acc	Mr Jatni Mahajan Executive Director	A 180-181, Okhla Industrial Area Phase-1, New Delhi 110 020 India Tel: +91-11-26818971-73	jmitra@mitra4u.com	www.mitra4u.com
	Advantage PAN Malaria Card	P	pLDH	Cassette	25/50/100 tests/pack Individual	No	24	4-30 C, Acc				
	Advantage Malaria Sixx	F, V	HRP2/pLDH	Dipstick	30/50/100 tests/pack Individual	Yes	24	4-30 C, Acc				
	Advantage MAL Sixx	P, F	pLDH	Dipstick	20/50/100 tests/pack Individual	Yes	24	4-30 C, Acc				
Merlin Labs, Inc.	Malaria Combo	F, V, O, M, P	HRP2 Plasmodium aldolase	Dipstick, Cassette	Individual	Yes	25	Acc	Mr Robert Hudak, Owner Ms Nicole Burns Management Representative	8082 Carle Del Cielo Carlsbad, CA 92011 U.S.A. Tel: 760-804-1782 Fax: 760-804-3885	rhudak@merlinlabs.com	www.merlinlabs.com
	Malaria P.I.	PI	HRP2	Dipstick, Cassette	Individual	Yes	25	Acc				
	Malaria Fever	F, V, O, M, P	Plasmodium aldolase	Dipstick, Cassette	Individual	Yes	25	Acc				

LIST OF KNOWN COMMERCIALY-AVAILABLE ANTIGEN-DETECTING MALARIA RDTs
Information for national public health services and UN Agencies wishing to procure RDTs
 (This list is not based on any criteria of product quality)

See notes on purchase and use of RDTs to guide use of this list, <http://www.wpro.who.int/rdt/link4.cfm>

Disclaimer: Inclusion in this list does not in any way imply an endorsement, certification, warranty of fitness or recommendation by WHO of any company or product for any purpose, and does not imply preference over products of a similar nature that are not mentioned. WHO does not furthermore warrant that: (1) the list is complete and/or that the information provided is or will be always accurate, and/or that (2) the products listed are of acceptable quality, have obtained regulatory approval in any country, or that their use is otherwise in accordance with the national laws and regulations of any country, including but not limited to patent laws. The list is based on evidence of quality manufacturing, through provision of evidence of independent certification of compliance with ISO13485:2003 (or US FDA 21 CFR 820) and a stability testing protocol provided by the companies to WHO, and whilst requesting this evidence from manufacturers as a pre-requisite for listing of products. WHO can not guarantee the accuracy of this information. WHO disclaims any and all liability and responsibility for any injury, death, loss, damage, or other prejudice of any kind whatsoever that may arise as a result of, or in connection with the procurement, distribution and use of any product included on the list. WHO recommends that the products should only be used by individuals trained

Products	Product name	Plasmodium species targeted (F = P falciparum V = P vivax O = P ovale M = P malariae P = pan-specific)	Target Antigen	Format	Packaging	All materials can be included? (Lancets, swabs, wells if required)	Shelf Life (months)	Storage temperature and stability data on which Shelf Life is based (Acc = Accelerated or RT = Real Time)	Contact	Address & Contact Nos.	Email	URL
Omega Diagnostics Ltd.	Visitect Malaria	F	HRP2	Cassette	Individual	Yes	24	4-30 C, RT	Mr Andrew Shephard Managing Director	Omega House Hillfoots Business Village Aha, Clacks, FK12 8DD Scotland, United Kingdom Tel: +44 1259 763030 Fax: +44 1259 761853	odf@omegadiagnostics.co.uk andrews@omegadiagnostics.co.uk	www.omegadiagnostics.co.uk
	Visitect Malaria Combo	F, P	HRP2/pLDH	Cassette	Individual	Yes	24	4-30 C, RT				
Orchid Biomedical Systems	ParaCheck	F	HRP2	Dipstick	25 tests/box Individual	Yes	24	4-40 C, RT	N. Sriam Director/Partner	Plot Nos 8889, Phase IIC Verna Industrial Estate Verna Goa 403 722, India Tel: +91 832 278 3140 Fax: +91 832 278 3139/2139	orchid@tulipgroup.com orchid_goa@sancharnet.in	www.tulipgroup.com
	ParaCheck	F	HRP2	Cassette	25 tests/box	Yes	24	4-40 C, RT				
Premier Medical Corporation	First Response® Ag	F, V	pLDH	Cassette	30 tests/box	Yes	18	2-30 C, Acc	Mr Neil Mehta President and CEO	1304 Johnston Drive Watchung NJ 07069 Tel: +1-908-344-6266 Fax: +1-530-869-7966 Mobile: +1-732-586-4422	nleah@premiermedcorp.com www.premiermedcorp.com	www.premiermedcorp.com
	First Response Ag	F	HRP2	Cassette	30 tests/box	Yes	18	2-30 C, Acc				
PT INDEC Diagnostics	Plasmodic Malaria PI	F, P	HRP2/pLDH	Cassette	25 tests/box, Individual	Sample loop Clearing buffer	24	2-30 C, RT	Reni Pangastuti	Taman Modern A4-2 Jl. Raya Bekasi Km. 24 Jakarta (13860), Indonesia Tel: 62 21 460 8922 Fax: 62 21 460 9977	corp@indec-diagnostics.co.id reni_pangastuti@indec-	
	Plasmodic Malaria P9Pv	F, V	HRP2/pLDH	Cassette	25 tests/box, Individual	Sample loop Clearing buffer	18	2-30 C, RT				
	Plasmodic Malaria-3	F, P, V	HRP2/pLDH	Cassette	25 tests/box, Individual	Sample loop Clearing buffer	24	2-30 C, RT				
Span Diagnostics Ltd.	ParasHT-4 Dipstick Test	F	HRP2	Dipstick	10, 50 tests/box, Individual	Sterile lancets, Isopropyl swabs and blood collection capillary and test tubes	24	Below 30C, RT	Mr Khushroo Pastakia Manager-International Business	173/B, New Industrial Estate, Road #6C, Udina, Surabaya 394210, India Tel. no. +91 261 2277211 Fax no. +91 261 2278319	export@spandag.com	www.span.in
	ParasHT-4 Device Test	F	HRP2	Cassette	10, 50 tests/box, Individual	Sterile lancets, Isopropyl swabs and blood collection capillary	24	Below 30C, RT				
	ParasHT-Total Dipstick Test	F, P	HRP2/aldolase	Dipstick	10, 50 tests/box, Individual	Sterile lancets, Isopropyl swabs and blood collection capillary and	24	Below 30C, RT				
	ParasHT-Total Device Test	F, P	HRP2/aldolase	Cassette	10, 50 tests/box, Individual	Sterile lancets, Isopropyl swabs and blood collection capillary	24	Below 30C, RT				
Standard Diagnostics Inc.	SD Bioline Malaria Ag (Antigen Test)	F, P	pLDH	Cassette	25 tests/box, Individual	Yes	24	1-40 C, Acc and RT	Mr John H. Kim VP, Sales and Marketing Manager/Sales	156-68 Hagal-dong, Gheungju, Yonjin-si, Kyonggi-do, Korea 446-930 Tel. no. 82-31-899-9700/9732 Fax no. 82-31-899-9740	johnhkim@standardia.com joyle@standardia.com sales@standardia.com	www.standardia.com
	SD Bioline Malaria Ag PI	F	HRP2	Cassette	25 tests/box, Individual	Yes	24	1-40 C, Acc and RT				
	SD Bioline Malaria Ag P/PPan	F, P	pLDH/HRP2	Cassette	25 tests/box, Individual	Yes	24	1-40 C, Acc and RT				
Unimed International, Inc.	FirstSign Malaria	F	HRP2	Cassette	Individual or Bulk	Yes	24	1-40 C, RT	Dr Arun Jaiswal President	2950 Scott Boulevard Santa Clara California 95054 U.S.A Tel: (408) 962-9156 Fax: (408) 962-9356	unimedinc@aol.com info@unimed.net	www.unimed.net
	FirstSign Malaria	F	HRP2	Dipstick	Individual or Bulk	Yes	24	4-30 C, RT				
	FirstSign PanCheck	P	pLDH	Cassette	Individual or Bulk	Yes	24	4-30 C, RT				
	FirstSign PanCheck	P	pLDH	Dipstick	Individual or Bulk	Yes	24	4-30 C, RT				
	FirstSign Para View	F, P	pLDH/HRP2	Cassette	Individual or Bulk	Yes	24	4-30 C, RT				
	FirstSign Para View -2	F, V	pLDH/HRP2	Cassette	Individual or Bulk	Yes	24	4-30 C, RT				
	FirstSign Para View -3	F, V, P	pLDH/HRP2	Cassette	Individual or Bulk	Yes	24	4-30 C, RT				
Vision Biotech	Malaria Rapid PI	F	HRP2	Cassette	25 or 100 tests/box	Yes	24	4-28 C, RT	Nick Borain, Director	Vision Biotech, 5, Thomas Close, 7 Inyoni Street, Ndabeni, Cape Town, South Africa	nickb@visionbiotech.com	
	Malaria Rapid Combo	P, F	HRP2/Aldolase	Cassette	25 or 100 tests/box	Yes	24	4-28 C, RT				
	Malaria Rapid Dual	P, F	HRP2/pLDH	Cassette	25 or 100 tests/box	Yes	24	4-28 C, RT				
Zephyr Biomedical Systems	Parascreen Device	P, F	pLDH, HRP2	Cassette	1, 5, 10, 25, 50, 100 Individual tests/box	Yes	24	4-30 C, RT and Acc	N. Sriam Director/Partner	Plot Nos M46/47, Phase IIIB Verna Industrial Estate Verna Goa 403 722, India Tel: +91 832 288 7023, 288 7025 Fax: +91 832 288 7024/2139	zephyr@tulipgroup.com	www.tulipgroup.com
	Falovax Device	F, V	pLDH, HRP2	Cassette	1, 5, 10, 25, 50, 100 Individual tests/box	Yes	24	4-30 C, RT and Acc				
	Parimax 3 Device	P, V, F	pLDH, HRP2	Cassette	1, 5, 10, 25, 50, 100 Individual tests/box	Yes	24	4-30 C, RT and Acc				
	Parabank Device	P	pLDH	Cassette	1, 5, 10, 25, 50, 100 Individual tests/box	Yes	24	4-30 C, RT and Acc				
	Parabank Dipstick	P	pLDH	Dipstick	1, 5, 10, 25, 50, 100 Individual tests/box	Yes	24	4-30 C, RT and Acc				
	MelaScan Device	PI, P	HRP2, Pan specific, aldolase	Cassette	1, 5, 10, 25, 50, 100 Individual tests/box	Yes	24	4-30 C, Acc				